


# INDIBA

## DEEP BEAUTY

PROIONIC BODY CARE SYSTEM®



# Visible results from the first applications\*



Control pictures taken before and after Indiba Deep Beauty treatments. Original image with no retouching.

5 APPLICATIONS



5 APPLICATIONS




5 APPLICATIONS

\*The results are related to 3 basic factors: the number of sessions, the energy applied and the volume of the tissue.

# INDIBA

## DEEP BEAUTY

PROIONIC BODY CARE SYSTEM®



3 different models to fit your specific needs.  
Capacitive and resistive modes.

ER 40 Cap: 250 VA Res: 65 W

ER 42 Cap: 350 VA Res: 100 W

ER 45 Cap: 450 VA Res: 200 W

## Indiba conductive gels

Developed to ensure the optimal transfer to the body for the perfect application of Indiba Deep Beauty system.

Moisturizing Gel · Anti cellulite Gel · Eye contour Gel

.....  
Contact us for additional technical information.

**INDIBA UK** International House

Southampton International Business Park

George Curl Way

Southampton

SO18 2RZ

United Kingdom

Tel.: +44 (0) 0845 4680560

Fax: +44 (0) 0238 0302157

[www.indiba-care.co.uk](http://www.indiba-care.co.uk)

[enquires@indiba-care.co.uk](mailto:enquires@indiba-care.co.uk)

The images and texts of this document and the trademark Indiba Deep Beauty Proionic Body Care System are protected and registered. No reproduction allowed in whole or in part.

MKT-2010-009 V01-2010-02-18 EN